

## Subsections

- [Labels](#)
  - [Example of a built-in analysis tool: the RMSD Trajectory Tool](#)
  - [Example of an analysis script](#)
- 

# Data Analysis in VMD

VMD is a powerful tool for analysis of structures and trajectories. Numerous tools for analysis are available under the VMD Main menu item *Extensions* → *Analysis*. In addition to these built-in tools, VMD users often use custom-written scripts to analyze desired properties of the simulated systems. VMD Tcl scripting capabilities are very extensive, and provide boundless opportunities for analysis. In this following section, we will learn how to use built-in VMD features for standard analysis, as well as consider a simple example of scripting.

## Labels

Labels can be placed in VMD to get information on a particular selection, to be used during visualization and quantitative analysis. Labels are selected with the mouse and can be accessed in *Graphics* → *Labels* menu. We will cover labels that can be placed on atoms and bonds, although angle and dihedral labeling are also available. In this context, labels for ``bonds" or ``angles" actually mean distances between two atoms or angles between three atoms; the atoms do not have to be physically connected by bonds in the molecule.

1

Start a new VMD session. Load the ubiquitin trajectory into VMD (using the files `ubiquitin.psf` and `pulling.dcd`). For graphical representation, display protein only, using *New Cartoon* for drawing method and *Structure* for coloring method. If you need help, check Section [2.1](#), steps 1-5.

2

Choose the *Mouse* → *Labels* → *Atoms* menu item from the VMD Main menu. The mouse is now set to the mode for displaying atom labels. You can click on any atom on your molecule and a label will be placed for this atom. Clicking again on it will erase the label.

3

We will now try the same for bonds. Choose the *Mouse* → *Label* → *Bonds* menu item from the VMD Main menu. This selects the ``Display Label for Bond" mode.


We will consider the distance between the  $\alpha$  carbon of Lysine 48 and of the C terminus. In the pulling simulation, the former is kept fixed, and the latter is pulled at a constant force of 500 pN. In reality, polyubiquitin chains can be linked by a connection between the C terminus of one ubiquitin molecule and the Lysine 48 of the next. The simulation then mimics the effect of pulling on the C terminus with this kind of linkage.

4

We will make a VDW representation for the  $\alpha$  carbons of Lysine 48 and of the C terminus. To find out the index of these atoms, make a selection including these two atoms, by typing in the Tk Console window (to


open the Tk Console window, select *Extensions* → *Tk Console* in the VMD Main menu.):  
set sel [atomselect top "resid 48 76 and name CA"].

- 5
- Get the indices by typing the following line in the Tk Console window:  
\$sel get index  
This command should give the indices 770 1242.




- 6
- In the Graphical Representations window, create a representation for the selection index 770 1242, with *VDW* as drawing method.

- 7
- Now that you can see the two  $\alpha$  carbons, choose the *Mouse* → *Label* → *Bonds* menu item from the VMD Main menu. Click on each atom one after the other. You should get a line connecting the two atoms (Fig. 37). The number appearing next to the line is the distance between the two atoms in Ångstroms. Note that the appearance of the line (its color), as well as appearance of essentially all other objects in VMD, can be changed in *Graphics* → *Colors* in the VMD Main menu.



- 8
- The value of the distance displayed corresponds to the current frame. Try playing the trajectory - you will see that the label is modified automatically as the distance between the atoms changes.



**Figure 37:** Labels in VMD. Label control is available in the *Graphics* → *Labels* menu item, using which, one can, e.g., plot the labeled distance as a function of time.

9

The labels can be used not only for display, but also for obtaining quantitative information. In VMD Main menu, select *Graphics* → *Labels*. On the top left-hand side of the window, there is a pull-down menu where you can choose the type of label (*Atoms*, *Bonds*, *Angles*, *Dihedrals*). For now, keep it in *Atoms*. You can see the list of atoms for which you have made a label.

10

Click on one of the atoms. You can see all the information of the atom displayed on the bottom half of the Labels window. This information is useful to make selections; it corresponds to the current frame, and is updated as the frame is changed.

11

You can also delete, hide, or show the atom label by clicking on the corresponding button on the top of the Labels window.

12

Now, in the Labels window, choose the label type *Bonds*, and select the ``bond" (distance) you labeled (Fig. 37). The information given corresponds to only the first atom in the bond, but the number in the *Value* field corresponds to the length of the bond in Ångströms.

13


Click on the *Graph* tab. Select the bond you labeled between atoms 770 and 1242. Click on the *Graph* button. This will create a plot of the distance between these two atoms over time (Fig. 37). You can also save this data to a file by clicking on the *Save* button, and then use an external plotting program to visualize the data.

14

Quit VMD.

## Example of a built-in analysis tool: the RMSD Trajectory Tool

The built-in analysis tools in VMD are available under the menu item *Extensions* → *Analysis*. These tools each features a GUI window that allow one to enter parameters and customize the quantities analyzed. In addition, all tools can be invoked in a scripting mode, using the TkConsole window. We will learn how to work with one of the most frequently used tools, the *RMSD Trajectory Tool*.



In this example, we will analyze RMSD for two trajectories for the same system, `ubiquitin.psf`. One of them is the already familiar pulling trajectory, `pulling.dcd`, and the other is the trajectory of a simulation in which no force was applied to the protein, `equilibration.dcd`.

1


Start a new VMD session. Load the ubiquitin equilibration trajectory into VMD (using the files `ubiquitin.psf` and `equilibration.dcd`).

2

Choose *Extensions* → *Analysis* → *RMSD Trajectory Tool* in the VMD Main window (Fig. [38](#)). The RMSD Trajectory Tool window will show up.

3

In the RMSD Trajectory Tool window, you can see many customizations can be made. For the default values, the molecule to be analyzed is `ubiquitin.psf` (the only one loaded). The selection for which RMSD will be computed is all of the `protein` atoms, excluding hydrogens (since the `noh` checkbox is on). The RMSD will be calculated for each frame (time  $t_1$  in Eq. [1](#)) with the reference to frame 0 (time  $t_2$  in Eq. [1](#)). Make sure the *Plot* checkbox is selected.



**Figure 38:** RMSD Trajectory Tool. The RMSD is plotted for the equilibration of ubiquitin.

4

Click the *Align* button. This will align each frame of the trajectory with respect to the reference frame (in this case, frame 0) to minimize RMSD, by applying only rigid-body translations and rotations. This step is not necessary, but is desirable in most cases, because we are interested only in RMSD that arises from the fluctuations of the structure and not from the displacements and rotations of the molecule as a whole. The result of the alignment can be seen in the OpenGL display.

5

Now, click the *RMSD* button in the RMSD Trajectory Tool window. The protein RMSD (in Ångstroms) vs. frame number is displayed in a plot (Fig. [38](#)).



**RMSD plot.** Over the several initial frames,  $\text{RMSD} = 0$  because positions of the protein atoms are fixed during that time in the simulation, to allow water molecules around the protein to adjust to the protein surface. After that, the protein is released, and RMSD grows quickly up to around  $1.5 \text{ \AA}$ . At that point, RMSD levels off and remains at  $\sim 1.5 \text{ \AA}$  further on. This is a typical behavior for MD simulations. Leveling of the RMSD means that the protein has relaxed from its initial crystal structure (which is affected by crystal packing and usually misses some atoms, e.g., hydrogens) to a more stable one. Production MD simulations are usually preceded by such equilibration runs, where the protein is allowed to relax; the process is monitored by checking RMSD vs. time, and equilibration is assumed to be sufficient when RMSD levels off. The RMSD of  $1.5 \text{ \AA}$  is quite a good value, totally acceptable for most protein simulations. Usually, the deviations from the crystal structure in a simulation are due to the thermal motion and to the relaxation process mentioned; imperfections of the simulation force-fields can contribute as well.

6

We will now work with the other trajectory. Load the pulling trajectory into VMD using the files `ubiquitin.psf` and `pulling.dcd`. Make sure you load `ubiquitin.psf` as a new molecule. You can change the names of the molecules by double-clicking on them in the VMD Main menu (see Section [4.1.2](#)).

7


In the RMSD Trajectory Tool window, hit the button *Add all* to update the list of molecules.

8

Click the *Align* button, and then click *RMSD*. The new graph (Fig. [39](#)) displays two RMSD plots vs. time, one for the equilibration trajectory, and the other for the pulling trajectory. The pulling RMSD does not level off and is much higher than the equilibration RMSD, since the protein is stretched in the simulation.

9

Quit VMD.



**Figure 39:** RMSD vs. time for the equilibration (blue) and pulling (red) trajectories of ubiquitin.

## Example of an analysis script

In many cases, one requires special types of trajectory analysis that are tailored for certain needs. The Tcl scripting in VMD provides opportunities for such custom tasks. Users commonly write their own scripts to analyze the features of interest, and a very extensive library of VMD scripts, contributed by many users, is available online, at [http://www.ks.uiuc.edu/Research/vmd/script\\_library/](http://www.ks.uiuc.edu/Research/vmd/script_library/).

We will investigate a very simple exemplary script, `distance.tcl`, which computes the distance between two atom selections vs. time and also distribution of the distance.

1

Start a new VMD session. Load the ubiquitin equilibration trajectory (files `ubiquitin.psf` and `equilibration.dcd`).

2

Open the TkConsole window by selecting *Extensions* → *Tk Console* in the VMD Main menu.

3

In the the TkConsole window, load the script into VMD by typing: `source distance.tcl` (make sure that the file `distance.tcl` is in the current folder). This will load the procedure defined in `distance.tcl` into the VMD.

4

One can now invoke the procedure by typing `distance` in the the TkConsole window. In fact, the correct usage is

```
distance seltext1 seltext2 N_d f_r_out f_d_out
```

where `seltext1` and `seltext2` are the selection texts for the groups of atoms between which the distance is measured, `N_d` is the number of bins for the distribution, and `f_r_out` and `f_d_out` are the file names to where the output distance vs. time and distance distribution will be written.

## 5

Open the script file `distance.tcl` with a text editor. You can see that the script does the following:

- Choose atom selections
 

```
set sel1 [atomselect top "$seltext1"]
set sel2 [atomselect top "$seltext2"]
```
- Get the number of frames in the trajectory and assign this value to the variable `nf`

```
set nf [molinfo top get numframes]
```
- Open file specified by the variable `f_r_out`

```
set outfile [open $f_r_out w]
```
- Loop over all frames
 

```
for {set i 0} {$i < $nf} {incr i} {
```
- Write out the frame number and update the selections to the current frame
 

```
puts "frame $i of $nf"
$sel1 frame $i
$sel2 frame $i
```
- Find the center of mass for each selection (`com1` and `com2` are position vectors)
 

```
set com1 [measure center $sel1 weight mass]
set com2 [measure center $sel2 weight mass]
```
- At each frame `i`, find the distance by subtracting one vector from the other (command `vecsub`) and computing the length of the resulting vector (command `veclength`), assign that value to an array element `simdata($i.r)`, and print a frame-distance entry to a file
 

```
set simdata($i.r) [veclength [vecsub $com1 $com2]]
puts $outfile "$i $simdata($i.r)"
}
```
- Close the file
 

```
close $outfile
```
- The second part of the script is for obtaining the distance distribution. It starts from finding the maximum and minimum values of the distance
 

```
set r_min $simdata(0.r)
set r_max $simdata(0.r)
for {set i 0} {$i < $nf} {incr i} {
  set r_tmp $simdata($i.r)
  if {$r_tmp < $r_min} {set r_min $r_tmp}
  if {$r_tmp > $r_max} {set r_max $r_tmp}
}
```
- The step over the range of distances is chosen based on the number of bins `N_d` defined in the beginning, and all values for the elements of the distribution array are set to zero
 

```
set dr [expr ($r_max - $r_min) / ($N_d - 1)]
for {set k 0} {$k < $N_d} {incr k} {
  set distribution($k) 0
}
```
- The distribution is obtained by adding 1 (`incr ...`) to an array element every time the distance is within the respective bin
 

```
for {set i 0} {$i < $nf} {incr i} {
```



```

set k [expr int(($simdata($i.r) - $r_min) / $dr)]
incr distribution($k)
}


```

- Write out the file with the distribution

```

set outfile [open $f_d_out w]
for {set k 0} {$k < $N_d} {incr k} {
puts $outfile "[expr $r_min + $k * $dr] $distribution($k)"
}
close $outfile

```



**Figure 40:** Distance between a residue and the center of ubiquitin. The distances analyzed are those for residue 76 (black) and residue 10 (green).

Now run the script by typing in the TkConsole window

```
distance "protein" "protein and resid 76" 10 res76-r.dat res76-d.dat
```

This will compute the distance between the center of the protein and center of the terminal residue 76, and write the distance vs. time and its distribution to files `res76-r.dat` and `res76-d.dat`.

7

Repeat the same for the protein's residue 10 by typing in the TkConsole window

```
distance "protein" "protein and resid 10" 10 res10-r.dat res10-d.dat
```

The data in files produced by the script `distance.tcl` are in two-column format. Compare the outputs for residue 76 and 10 using your favorite external plotting program (Fig. [40](#)).

Residue 76 is at the protein's C-terminus, which is extended towards the solvent and is quite flexible, while residue 10 is at the surface of the globular part of ubiquitin. The difference in their dynamics with respect to the rest of the protein is immediately obvious when our newly obtained data are plotted (Fig. [40](#)): the distance of residue 76 from the protein's center is substantially greater than that of residue 10, and the distribution of the distance is noticeably wider due to the flexibility of the C-terminus. This is just a simple example of scripting for the analysis of a trajectory. Similar, but usually much more complex, customized scripts are routinely employed by VMD users to perform many kinds of analysis.

8

Quit VMD.

---

[Next](#) [Up](#) [Previous](#)

**Next:** [Running VMD on Supercomputers](#) **Up:** [VMD Tutorial](#) **Previous:** [Comparing Structures and Sequences](#)  
*vmd@ks.uiuc.edu*